


Product Name: Recombinant Human AIRE Protein
Catalog #: PEH10470



Summary

Name	AIRE_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43918
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ALVPPRLPT KRKASEEARA AAPAALTPRG TASPGSQLKA KPPKKPESSA EQQRLPLGNG IQTMSASVQR AVAMSSGDVP GARGAVEGIL IQQVFESGGS KKCIQVGGEF YTPSKFEDSG SGKNKARSSS GPKPLVRAKG AQGAAPGGGE ARLGQQGSVP APLALPSDPQ LHQKNEDECA

image



Background

Alternative Names Autoimmune regulator (Autoimmune polyendocrinopathy candidiasis ectodermal dystrophy protein) (APECED protein)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human AIRE Protein
Catalog #: PEH10470

