


Product Name: Recombinant Human GLI3 Protein
Catalog #: PEH10472



Summary

Name	GLI3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P10071
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QQEQPEGTTL VKEEGDKDES KQEPEVIYET NCHWEGCARE FDTQEQLVHH INNDHIHGEK KEFVCRWLDC SREQKPFKAQ

image



Background

Alternative Names	Transcriptional activator GLI3 (GLI3 form of 190 kDa) (GLI3-190) (GLI3 full-length protein) (GLI3FL) [Cleaved into: Transcriptional repressor GLI3R (GLI3 C-terminally truncated form) (GLI3 form of 83 kDa) (GLI3-83)]
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use