


Product Name: Recombinant Human PYCARD Protein
Catalog #: PEH10480



Summary

Name	ASC_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9ULZ3
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MGRARDAILD ALENLTAEEL KFKLKLLSV PLREGYGRIP RGALLSMDAL DLTDKLVSYF LETYGAELTA NVLRDMGLQE MAGQLQAATH QGSGAAPAGI QAPPQSAAKP GLHFIDQHRA ALIARVTNVE WLLDALYGKV LTDEQYQAVR AEPTNPSKMR KLFSTPAWN WTCKDLLLQA LRESQSYLVE DLERS

image



Background

Alternative Names	Apoptosis-associated speck-like protein containing a CARD (hASC) (Caspase recruitment domain-containing protein 5) (PYD and CARD domain-containing protein) (Target of methylation-induced silencing 1)
--------------------------	---

Note

Product Name: Recombinant Human PYCARD Protein
Catalog #: PEH10480



For Research Use Only , Not for Diagnostic Use