


Product Name: Recombinant Human E2F1 Protein
Catalog #: PEH10483



Summary

Name	E2F1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q01094
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VDLNWAAEVL KVKRRRIYDI TNVLEGIQLI AKKSKNHIQW LGSHTTVGVG GRLEGLTQDL RQLQESEQQL DHLMNICTTQ LRLLEDTDS QRLAYVTCQD LRSIADPAEQ MVMVIKAPPE TQLQAVDSSE NFQISLKSQK GPIDVFLCPE ETVGGISPGK TPSQEVTSSEE ENRATDSATI VSPPPSSPPS

image



Background

Alternative Names	Transcription factor E2F1 (E2F-1) (PBR3) (Retinoblastoma-associated protein 1) (RBAP-1) (Retinoblastoma-binding protein 3) (RBBP-3) (pRB-binding protein E2F-1)
--------------------------	---

Note

Product Name: Recombinant Human E2F1 Protein
Catalog #: PEH10483



For Research Use Only , Not for Diagnostic Use