



Product Name: Recombinant Human AMPH Protein
Catalog #: PEH10486



Summary

Name	AMPH_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P49418
Host	E.coli
Species	Human
Molecular Mass	43 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KNVQKRLNRA QEKVLQKLGK ADETKDEQFE EYVQNFKRQE AEGTRLQREL RGLAAIKGM QEASMKLTES LHEVYEPDWY GREDVKMVGE KCDVLWEDFH QKLVDGSLLT LDTYLGQFPD IKNRIAKRSR KLVDYDSARH HLEALQSSKR KDESRISKAE EEFAQKQKVF EEFNVDLQEE LPSLWSRRVG FYVNTFKNVS SLEAKFHKEI AVLCHKLYEV MTKLGDQHAD KAFTIQGAPS DSGPLRIAKT PSPPEEPSPL PSPTASPNHT LAPASPAPAR PRSPSQTRKG PPVPPLPKVT PTKELQQENI ISFFEDNFVP EISVTTPSQN EVPEVKKEET LLDLDFDPFK PEVTPAGSAG

image



Background

Alternative Names	Amphiphysin
--------------------------	-------------

Product Name: Recombinant Human AMPH Protein
Catalog #: PEH10486



Note

For Research Use Only , Not for Diagnostic Use