


Summary

Name	AMPH_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P49418
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EVPEVKKEET LLDLDFDPFK PEVTPAGSAG VTHSPMSQTL PWDLWTTSTD LVQPASGGSF NGFTQPQDTS LFTMQTDQSM ICNLAESEQA PPTEPKAEEP LAAVTPAVGL DLGMDTRAEE PVEEAVIIPG ADADAAVGTL VSAAE GAPGE EAEA EKATVP AGEGVSLEEA KIGTETTEGA ESAQPEAEEL EATVPQEKVI PSVVI EPASN HEEEGENEIT IGAEPKETTE DAAPPGPTSE TPELATEQKP IQDPQPTPSA PAMGAADQLA SAREASQELP PGFLYKVELT HDFEAANSDE LTLQRGDVVL VVPSDSEADQ DAGWL VGVKE SDWLQYRDLA TYKGLFPENF

image



Background

Alternative Names Amphiphysin

Product Name: Recombinant Human AMPH Protein
Catalog #: PEH10487



Note

For Research Use Only , Not for Diagnostic Use