


Product Name: Recombinant Human AOX1 Protein
Catalog #: PEH10490



Summary

Name	AOXA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q06278
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QRTRVFGSER MMWFSPVTLK ELLEFKFKYP QAPVIMGNTS VGPEVKFKGV FHPVIISPDR IEELSVVNHA YNGLTLGAGL SLAQVKDILA DVVQKLPEEK TQMYHALLKH LGTLAGSQIR NMASLGGHII SRHPDSLNP ILAVGNCTLN LLSKEGKRQI PLNEQFLSKC PNADLKPQEI LSVNIPYSR KWEFVSAFRQ

image



Background

Alternative Names Aldehyde oxidase (EC 1.2.3.1) (Aldehyde oxidase 1) (Azaheterocycle hydroxylase) (EC 1.17.3.-)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human AOX1 Protein
Catalog #: PEH10490

