


Product Name: Recombinant Human EEA1 Protein
Catalog #: PEH10498



Summary

Name	EEA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15075
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KTAKATLEQD SAKKEQQLQE RCKALQDIQK EKSLKEKELV NEKSKLAEIE EIKCRQEKEI TKLNEELKSH KLESIKEITN LKDAKQLLIQ QKLELQGKAD SLKAAVEQEK RNQQILKDQV KKEEEELKKE FIEKEAKLHS EIKEKEVGMK KHEENEAKLT MQITALNENL GTVKKKEWQSS

image



Background

Alternative Names	Early endosome antigen 1 (Endosome-associated protein p162) (Zinc finger FYVE domain-containing protein 2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human EEA1 Protein
Catalog #: PEH10498

