


Product Name: Recombinant Human TRADD Protein
Catalog #: PEH10499



Summary

Name	TRADD_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15628
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	WVGSAFLVE SSDKVVLSD AYAHPQKVA VYRALQAALA ESGGSPDVLQ MLKIHRSDPQ LIVQLRFCGR QPCGRFLRAY REGALRAALQ RSLAAALAQH SVPLQLELRA GAERLDALLA DEERCLSCIL AQQPDRLRDE ELAELEDALR NLKCGSGARG

image



Background

Alternative Names	Tumor necrosis factor receptor type 1-associated DEATH domain protein (TNFR1-associated DEATH domain protein) (TNFRSF1A-associated via death domain)
--------------------------	--

Note

Product Name: Recombinant Human TRADD Protein
Catalog #: PEH10499



For Research Use Only , Not for Diagnostic Use