


**Product Name: Recombinant Human BIRC3 Protein**  
**Catalog #: PEH10506**



## Summary

<b>Name</b>	BIRC3_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q13489
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	14 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	YLIRIKQEF IRQVQASYPH LLEQLLSTSD SPGDENAESS IIHFEPGEDH SEDAIMMNTP VINAAVEMGF SRSLVKQTVQ RKILATGENY RLVNDLVLDL LNAEDEIREE

## image



## Background

<b>Alternative Names</b>	Baculoviral IAP repeat-containing protein 3 (EC 2.3.2.27) (Apoptosis inhibitor 2) (API2) (Cellular inhibitor of apoptosis 2) (C-IAP2) (IAP homolog C) (Inhibitor of apoptosis protein 1) (hIAP-1) (hIAP1) (RING finger protein 49) (RING-type E3 ubiquitin transferase BIRC3) (TNFR2-TRAF-signaling complex protein 1)
--------------------------	--

## Note

**Product Name: Recombinant Human BIRC3 Protein**  
**Catalog #: PEH10506**



---

For Research Use Only , Not for Diagnostic Use