

Product Name: Recombinant Human KHDRBS1 Protein
Catalog #: PEH10512



Summary

Name	KHDR1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q07666
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ATGPDATVGG PAPTLLPPS ATASVKMEPE NKYLPELMAE KDSLDPSTFH AMQLLTAEIE KIQKGDSSKD DEENYLDLFS HKNMCLKERV LIPVKQYPKF NFVGKILGPQ GNTIKRLQEE TGAKISVLGK GSMRDKAKEE ELRKGGDPKY

image



Background

Alternative Names	KH domain-containing, RNA-binding, signal transduction-associated protein 1 (GAP-associated tyrosine phosphoprotein p62) (Src-associated in mitosis 68 kDa protein) (Sam68) (p21 Ras GTPase-activating protein-associated p62) (p68)
--------------------------	--

Note

Product Name: Recombinant Human KHDRBS1 Protein
Catalog #: PEH10512



For Research Use Only , Not for Diagnostic Use