


Product Name: Recombinant Human EPS15 Protein
Catalog #: PEH10513



Summary

Name	EPS15_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P42566
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GKAQLEPLQQ HLQDSQQEIS SMQMKLMEMK DLENHNSQLN WCSSPHSILV NGATDYCSLS TSSSETANLN EHVEGQSNLE SEPIHQESPA RSSPELLPSG VTDENEVTTA VTEKVCSELD NNRHSKEEDP FNVDSSSLTG

image



Background

Alternative Names Epidermal growth factor receptor substrate 15 (Protein Eps15) (Protein AF-1p)

Note

For Research Use Only , Not for Diagnostic Use