


Product Name: Recombinant Human MSI1 Protein
Catalog #: PEH10521



Summary

Name	MSI1H_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43347
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ASPDSPHDPC KMFIGGLSWQ TTQEGLREYF GQFGEVKECL VMRDPLTKRS RGFGFVTFMD QAGVDKVLAQ SRHELDSTI DPKVAFPRRA QPKMVTRTKK IFVGGLSVNT TVEDVKQYFE QFGKVDDAML MFDKTTNRHR GFGFVTFESE DIVEKVCEIH FHEINNMVE CKKAQPKEVM SPTGSARGRS

image



Background

Alternative Names RNA-binding protein Musashi homolog 1 (Musashi-1)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MSI1 Protein
Catalog #: PEH10521

