


Product Name: Recombinant Human IGF2BP2 Protein
Catalog #: PEH10527



Summary

Name	IF2B2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y6M1
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FFNPKEEVKL EAHIRVPSST AGRVIGKGGK TVNELQNLTS AEVIVPRDQT PDNEEVIVR IIGHFFASQT AQRKIREIVQ QVKQQEQKYP

image



Background

Alternative Names	Insulin-like growth factor 2 mRNA-binding protein 2 (IGF2 mRNA-binding protein 2) (IMP-2) (Hepatocellular carcinoma autoantigen p62) (IGF-II mRNA-binding protein 2) (VICKZ family member 2)
--------------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Note

For Research Use Only , Not for Diagnostic Use