


Product Name: Recombinant Human DIDO1 Protein
Catalog #: PEH10530



Summary

Name	DIDO1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BTC0
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VDDLPNRMCA DVRRNSVERP AEPVAGAATP SLVEQQKMLE ELNKQIEEQK RQLEEQQEAL RQQRAAVGVS MAHFSVSDAL MSPPPKSSLP KAELFQQEQQ SADKPASLPP ASQASNHRDP RQARRLATET

image



Background

Alternative Names	Death-inducer obliterator 1 (DIO-1) (hDido1) (Death-associated transcription factor 1) (DATF-1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use