


Product Name: Recombinant Human DOCK1 Protein
Catalog #: PEH10532



Summary

Name	DOCK1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14185
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SMVRSFTMPS SSRPLSVASV SSLSSDSTPS RPGSDGFALE PLLPKKMHSR SQDKLDKDDL EKEKKDKKKE KRNSKHQEIF EKEFKPTDIS LQQSEAVILS ETISPLRPQR PKSQVMNVIG SERRFSVSPS SPSSQQTPPP VTPRAKLSFS MQSSLELNGM TGADVADVPP PLPLKGSVAD YGNLMENQDL LGSPTPPPPP PHQRHLPPPL PSKTPPPPPP KTTRKQASVD

image



Background

Alternative Names	Dedicator of cytokinesis protein 1 (180 kDa protein downstream of CRK) (DOCK180)
--------------------------	--

Note

Product Name: Recombinant Human DOCK1 Protein
Catalog #: PEH10532



For Research Use Only , Not for Diagnostic Use