


Product Name: Recombinant Human SIRT3 Protein
Catalog #: PEH10542



Summary

Name	SIR3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9NTG7
Host	E.coli
Species	Human
Molecular Mass	44 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GKLSLQDVAE LIRARACQRV VVMVGAGIST PSGIPDFRSP GSGLYSNLQQ YDLPYPEAIF ELPFFFHNP K PFFTLAKELY PGNYKPNVTH YFLRLLHDKG LLLRLYTQNI DGLERVSGIP ASKLVEAHGT FASATCTVCQ RPFPGEDIRA DVMADRVPRC PVCTGVVKPD IVFFGEPLPQ RFLHVVDFP MADLLLILGT SLEVFPFASL TEAVRSSVPR LLINRDLVGP LAWHPRSRDV AQLGDVVHGV ESLVELLGWT EEMRDLVQRE

image



Background

Alternative Names	NAD-dependent protein deacetylase sirtuin-3, mitochondrial (hSIRT3) (EC 2.3.1.286) (Regulatory protein SIR2 homolog 3) (SIR2-like protein 3)
--------------------------	--

Product Name: Recombinant Human SIRT3 Protein
Catalog #: PEH10542



Note

For Research Use Only , Not for Diagnostic Use