


Product Name: Recombinant Human TGFB11 Protein
Catalog #: PEH10545



Summary

Name	TGFI1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43294
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ASGEQKEDQS EDKKRPSLPS SPSGLPKAS ATSATLELDR LMASLSDFRV QNHLPASGPT QPPVVSSTNE GSPSPPEPTG KGSLDTMLGL LQSDLSRRGV

image



Background

Alternative Names Transforming growth factor beta-1-induced transcript 1 protein (Androgen receptor coactivator 55 kDa protein) (Androgen receptor-associated protein of 55 kDa) (Hydrogen peroxide-inducible clone 5 protein) (Hic-5)

Note

For Research Use Only , Not for Diagnostic Use