


Product Name: Recombinant Human ARRB1 Protein
Catalog #: PEH10546



Summary

Name	ARRB1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	P49407
Host	E.coli
Species	Human
Molecular Mass	42 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RNSVRLVIRK VQYAPERPGP QPTAETTRQF LMSDKPLHLE ASLDKEIYYH GEPISVNVHV TNNTNKT VKK IKISVRQYAD ICLFNTAQYK CPVAMEEADD TVAPSSTFCK VYTLTPFLAN NREKRG LALD GK LKHEDTNL ASSTLLREGA NREILGIIVS YKVKVKLVVS RG LLLGDLAS SDVAVELPFT LMHPKPKEEP PHREVPENET PVDTNLIELD TNDDDIVFED FARQLKGMK DDKEEEEDGT

image



Background

Alternative Names Beta-arrestin-1 (Arrestin beta-1) (Non-visual arrestin-2)

Note

Product Name: Recombinant Human ARRB1 Protein
Catalog #: PEH10546



For Research Use Only , Not for Diagnostic Use