

**Product Name: Recombinant Human SACS Protein**  
**Catalog #: PEH10548**



## Summary

|                              |  |
|------------------------------|--|
| <b>Name</b>                  | SACS_HUMAN   |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE  |
| <b>Endotoxin level</b>       | Please contact with the lab for this information   |
| <b>Tag</b>                   | N-MBP & C-His  |
| <b>Accession #</b>           | Q9NZJ4   |
| <b>Host</b>                  | E.coli   |
| <b>Species</b>               | Human  |
| <b>Molecular Mass</b>        | 108 kDa  |
| <b>Buffer</b>                | Phosphate buffered saline  |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.   |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.   |
| <b>Sequence</b>              | HHTAKLPADL RLSISVIDSS DEATIRLANM LKIEQLKTTS CLKLVLKDIE<br>NAFYSHHEVT QLMLWVLENL SSLKNENPNV LEWLTPLKFI QISQEQMVSA<br>GELFDPDIEV LKDLFCNEEG TYFPPSVFTS PDILHSLRQI GLKNEASLKE<br>KDVVQVAKKI EALQVGACPD QDVLLKKAKT LLLVLNKNHT LLQSSEGKMT<br>LKKIKWVPAC KERPPNYPGS LVWKGDLNCL CAPPDMCDVG HAILIGSSLP<br>LVESIHVNLE KALGIFTKPS LSAVLKHFKI VVDWYSSKTF SDEDYYQFQH<br>ILLEIYGFMH DHLNEGKDSF RALKFPVWWT GKKFCPLAQA VIKPIHDLDL<br>QPYLHNVPKT MAKFHQLFKV CGSIEELTSD HISMVIQKIY LKSDQDLSEQ<br>ESKQNLHLML NIIRWLYSNQ IPASPNTVPV IHHSKNPSKL IMKPIHECCY<br>CDIKVDDLND LLEDSEPII LVHEDIPMKT AEWLVKVPCLS TRLINPENMG<br>FEQSGQREPL TVRIKNILEE YPSVSDIFKE LLQNADDANA TECSFLIDMR<br>RNMDIRENLL |

## image



**SDS**

**Product Name: Recombinant Human SACS Protein**  
**Catalog #: PEH10548**



---

## Background

**Alternative Names**                      Sacsin (DnaJ homolog subfamily C member 29) (DNAJC29)

## Note

For Research Use Only , Not for Diagnostic Use