


Product Name: Recombinant Human GATA1 Protein
Catalog #: PEH10555



Summary

Name	GATA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P15976
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ARECVNCGAT ATPLWRRDRT GHYLCNACGL YHKMNGQNRP LIRPKKRLIV SKRAGTQCTN CQTTTTTLWR RNASGDPVCN ACGLYYKLHQ

image



Background

Alternative Names	Erythroid transcription factor (Eryf1) (GATA-binding factor 1) (GATA-1) (GF-1) (NF-E1 DNA-binding protein)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use