


Product Name: Recombinant Human TNRC6A Protein
Catalog #: PEH10558



Summary

Name	TNR6A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NDV7
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DVERNLSRDL VQEEQLMEE KKKKKDDKKK KEAAQKKATE QKIKVPEQIK PSVSQPQPAN SNGTSTATS TNNNAKRATA NNQPPQQQQQ QQQPQQQQPQ QQPQPQQQQ QPQQQPQALP RYPREVPPRF RHQEHKQLLK RGQHFPVIAA NLGSAVKVLN SQSESSALTN QQPQNNGEVQ NSKNQSDINH STSGSHYENS QRGPVSSTSD SSTNCKNAVV SDLSEKEAWP SAPGSDPELA SECMDADSAS SSESERNITI MASGNTGGEK DGLRNSTGLG SQNKFVVGSS SNNVGHGSST GPWGFSGHAI

image



Background

Alternative Names	Trinucleotide repeat-containing gene 6A protein (CAG repeat protein 26) (EMSY interactor protein) (GW182 autoantigen) (Protein GW1) (Glycine-
--------------------------	---

Product Name: Recombinant Human TNRC6A Protein
Catalog #: PEH10558



tryptophan protein of 182 kDa)

Note

For Research Use Only , Not for Diagnostic Use