


Product Name: Recombinant Human PREP Protein
Catalog #: PEH10559



Summary

Name	PPCE_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P48147
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RVFLDPNLS DDGTVALRGY AFSEEDGEYFA YGLSASGSDW VTIKFMKVDG AKELPDVLER VKFSCMAWTH DGKGMFYNSY PQQDGKSDGT ETSTNLHQKL YYHVLGTDQS EDILCAEFPD EPKWMGGAEL SDDGRYVLLS IREGCDPVNR LWYCDLQQES SGIAGILKWV KLIDNFEGEY DYVTNEGTVF

image



Background

Alternative Names Prolyl endopeptidase (PE) (EC 3.4.21.26) (Post-proline cleaving enzyme)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PREP Protein
Catalog #: PEH10559

