


Product Name: Recombinant Human SUMO3 Protein
Catalog #: PEH10560



Summary

Name	SUMO3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P55854
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSEEKPKGV KTENDHINLK VAGQDGSVVQ FKIKRHTPLS KLMKAYCERQ GLSMRQIRFR FDGQPINETD TPAQLEMEDE DTIDVFQQQT GG

image



Background

Alternative Names	Small ubiquitin-related modifier 3 (SUMO-3) (SMT3 homolog 1) (SUMO-2) (Ubiquitin-like protein SMT3A) (Smt3A)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use