


Product Name: Recombinant Human TOPBP1 Protein
Catalog #: PEH10564



Summary

Name	TOPBP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92547
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TNNKESAPSN GSGKNDKSGV LTQTLEMREN FQKQLQEIMS ATSIVKPQGQ RTSLSRSGCN SASSTPDSTR SARSGRSRVL EALRQSRQTV PDVNTEPSQN EQIIWDDPTA REERARLASN LQWPSCPTQY SELQVDIQNL EDSPFQKPLH DSEIAKQAVC DPGNIRVTEA

image



Background

Alternative Names	DNA topoisomerase 2-binding protein 1 (DNA topoisomerase II-beta-binding protein 1) (TopBP1) (DNA topoisomerase II-binding protein 1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TOPBP1 Protein
Catalog #: PEH10564

