


Product Name: Recombinant Human TNKS Protein
Catalog #: PEH10566



Summary

Name	TNKS1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O95271
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KHGADPTKKN RDGNTPLDLV KEGDTDIQDL LRGDAALLDA AKKGCLARVQ KLCTPENINC RDTQGRNSTP LHLAAGYNNL EVAEYLLEHG ADVNAQDKGG

image



Background

Alternative Names	Poly [ADP-ribose] polymerase tankyrase-1 (EC 2.4.2.30) (ADP-ribosyltransferase diphtheria toxin-like 5) (ARTD5) (Poly [ADP-ribose] polymerase 5A) (Protein poly-ADP-ribosyltransferase tankyrase-1) (EC 2.4.2.-) (TNKS-1) (TRF1-interacting ankyrin-related ADP-ribose polymerase) (Tankyrase I) (Tankyrase-1) (TANK1)
--------------------------	--

Note

Product Name: Recombinant Human TNKS Protein
Catalog #: PEH10566



For Research Use Only , Not for Diagnostic Use