


**Product Name: Recombinant Human HPGD Protein**  
**Catalog #: PEH10568**



## Summary

<b>Name</b>	PGDH_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P15428
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	30 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	TGAAQGIGRA FAEALLKGA KVALVDWNLE AGVQCKAALD EQFEPQKTLF IQCDVADQQQ LRDTFRKVVD HFGRLDILVN NAGVNNEKNW ECTLQINLVS VISGTYLGLD YMSKQNGGEG GIIINMSSLA GLMPVAQQPV YCASKHGIVG FTRSAALAN LMNSGVRLNA ICPGFVNTAI LESIEKEENM GQYIEYKDHI KDKIKYYGIL DPPLIANGLI TLIEDDALNG AIMKITTSGK IHFQDYDTTP

## image



## Background

<b>Alternative Names</b>	15-hydroxyprostaglandin dehydrogenase [NAD(+)] (15-PGDH) (EC 1.1.1.141) (Prostaglandin dehydrogenase 1) (Short chain dehydrogenase/reductase family 36C member 1)
--------------------------	---

**Product Name: Recombinant Human HPGD Protein**  
**Catalog #: PEH10568**

---



**Note**

For Research Use Only , Not for Diagnostic Use