


Product Name: Recombinant Human MAGI2 Protein
Catalog #: PEH10570



Summary

Name	MAGI2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86UL8
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EHDVRKPKEL SACGQKKQRL GEQRERSASP QRAARPRLEE APGGQGRPEA GRPASEARAP GLAAADAADA ARAGGKEAPR AAAGSELCRR

image



Background

Alternative Names	Membrane-associated guanylate kinase, WW and PDZ domain-containing protein 2 (Atrophin-1-interacting protein 1) (AIP-1) (Atrophin-1-interacting protein A) (Membrane-associated guanylate kinase inverted 2) (MAGI-2)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use