


Product Name: Recombinant Human CRYAA Protein
Catalog #: PEH10573



Summary

Name	CRYAA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P02489
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MDVTIQHPWF KRTLGPFPYS RLFDQFFGEG LFEYDLLPFL SSTISPYRQ SLFRTVLDSG ISEVRSDRDK FVIFLDVKHF SPEDLTVKVQ DDFVEIHGKH NERQDDHGYI SREFHRRYRL PSNVDQSALS CSLSADGMLT FCGPKIQTGL DATHAERAIP VSREEKPTSA PSS

image



Background

Alternative Names	Alpha-crystallin A chain (Heat shock protein beta-4) (HspB4) [Cleaved into: Alpha-crystallin A(1-172); Alpha-crystallin A(1-168); Alpha-crystallin A(1-162)]
--------------------------	--

Note

Product Name: Recombinant Human CRYAA Protein
Catalog #: PEH10573



For Research Use Only , Not for Diagnostic Use