


Product Name: Recombinant Human TRIM63 Protein
Catalog #: PEH10575



Summary

Name	TRI63_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q969Q1
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SRPLQKGSHP MCKEHEDEKI NIYCLTCEVP TCSMCKVFGI HKACEVAPLQ SVFQGQKTEL NNCISMLVAG NDRVQTIITQ LEDSRRTKE NSHQVKEELS QKFDTLYAIL DEKKSELLQR ITQEQEKKLS FIEALIQQYQ EQLDKSTKLV ETAIQSLDEP GGATFLLTAK QLIKSIVEAS KGCQLGKTEQ GFENMDFFTL DLEHIADALR AIDFGTDEEE EFIEEEDQE EEESTEGKEE

image



Background

Alternative Names	E3 ubiquitin-protein ligase TRIM63 (EC 2.3.2.27) (Iris RING finger protein) (Muscle-specific RING finger protein 1) (MuRF-1) (MuRF1) (RING finger protein 28) (RING-type E3 ubiquitin transferase TRIM63) (Striated muscle RING zinc finger protein) (Tripartite motif-containing protein 63)
--------------------------	---

Product Name: Recombinant Human TRIM63 Protein
Catalog #: PEH10575



Note

For Research Use Only , Not for Diagnostic Use