


Product Name: Recombinant Human RPS3 Protein
Catalog #: PEH10576



Summary

Name	RS3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P23396
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FVADGIFKAE LNEFLTRELA EDGYSGVEVR VTPTRTEIII LATRTQNVLG EKGRRIRELT AVVQKRFGFP EGSVELYAEK VATRGLCAIA QAESLRYKLL GGLAVRRACY GVLRFIMESG AKGCEVVVSG KLRGQRAKSM KFVDGLMIHS GDPVNYVVDV AVRHVLLRQG VLGKVKIML PWDPTGKIGP KKPLPDHVS VEPKDEILPT TPISQKGGK PEPPAMPQPV

image



Background

Alternative Names	40S ribosomal protein S3 (EC 4.2.99.18) (Small ribosomal subunit protein uS3)
--------------------------	---

Note

Product Name: Recombinant Human RPS3 Protein
Catalog #: PEH10576



For Research Use Only , Not for Diagnostic Use