


Product Name: Recombinant Human FRK Protein
Catalog #: PEH10586



Summary

Name	FRK_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P42685
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NPGALCSPQS QRHGHYFVAL FDYQARTAED LSFRAQDKLQ VLDTLHEGWW FARHLEKRRD GSSQQLQGYI PSNYVAEDRS LQAEPWFFGA IGRSDAEKQL LYSENKTGSF LIRESESQKG EFSLSVLDGA VVKHYRIKRL DEGGFFLTRR RIFSTLNEFV SHYTKTSDGL

image



Background

Alternative Names	Tyrosine-protein kinase FRK (EC 2.7.10.2) (FYN-related kinase) (Nuclear tyrosine protein kinase RAK) (Protein-tyrosine kinase 5)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human FRK Protein
Catalog #: PEH10586

