


**Product Name: Recombinant Human UBE2K Protein**  
**Catalog #: PEH10588**



## Summary

<b>Name</b>	UBE2K_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P61086
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	18 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	REFKEVLKSE ETSKNQIKVD LVDENFTELR GEIAGPPDTP YEGGRYQLEI KIPETYPFNP PKVRFITKIW HPNISSVTGA ICLDILKDQW AAAMTLRTVL LSLQALLAAA EPDDPQDAVV ANQYKQNPEN FKQTARLWAH

## image



## Background

<b>Alternative Names</b>	Ubiquitin-conjugating enzyme E2 K (EC 2.3.2.23) (E2 ubiquitin-conjugating enzyme K) (Huntingtin-interacting protein 2) (HIP-2) (Ubiquitin carrier protein) (Ubiquitin-conjugating enzyme E2-25 kDa) (Ubiquitin-conjugating enzyme E2(25K)) (Ubiquitin-conjugating enzyme E2-25K) (Ubiquitin-protein ligase)
--------------------------	---

**Product Name: Recombinant Human UBE2K Protein**  
**Catalog #: PEH10588**



---

**Note**

For Research Use Only , Not for Diagnostic Use