


Product Name: Recombinant Human CNDP2 Protein
Catalog #: PEH10589



Summary

Name	CNDP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96KP4
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IDENQDRYIK KLAKWVAIQS VSAWPEKRG E IRRMMEVAAA DVKQLGGSVE LVDIGKQKLP DGSEIPLPPI LLGRLGSDPQ KKTVCYGH L DVQPAALEDG WDSEPF TLVE RDGKLYGRGS TDDKGPVAGW INALEAYQKT GQEIPVNVRF CLEGMEESGS EGLDELIFAR KDTFFKDVDY VCISDNYWLG KKKPCITYGL RGICYFFIEV ECSNKDLHSG VYGGSVHEAM TDLILLMGSL VDKRGNILIP GINEAVAAVT EEEHKLYDDI DFDIEEFAKD

image



Background

Alternative Names	Cytosolic non-specific dipeptidase (EC 3.4.13.18) (CNDP dipeptidase 2) (Carnosine dipeptidase II) (Epididymis secretory protein Li 13) (Glutamate carboxypeptidase-like protein 1) (Peptidase A)
--------------------------	--

Product Name: Recombinant Human CNDP2 Protein
Catalog #: PEH10589



Note

For Research Use Only , Not for Diagnostic Use