


Product Name: Recombinant Human PIK3CG Protein
Catalog #: PEH10592



Summary

Name	PK3CG_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P48736
Host	E.coli
Species	Human
Molecular Mass	76 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DRVRAEMPNQ LRKQLEAIIA TDPLNPLTAE DKELLWHFRY ESLKHPKAYP KLFSSVKWGQ QEIVAKTYQL LARREVWDQS ALDVGLTMQL LDCNFSDENV RAIAVQKLES LEDDDVLHYL LQLVQAVKFE PYHDSALARF LLKRGLRNKR IGHFLFWFLR SEIAQSRHYQ QRFAVILEAY LRGCGTAMLH DFTQQVQVIE MLQKVTLDIK SLSAEKYDVS SQVISQLKQK LENLQNSQLP ESFRVPYDPG LKAGALAI EK CKVMASKKKP LWLEFKCADP

image



Background

Alternative Names	Phosphatidylinositol 4,5-bisphosphate 3-kinase catalytic subunit gamma isoform (PI3-kinase subunit gamma) (PI3K-gamma) (PI3Kgamma) (PtdIns-3-kinase subunit gamma) (EC 2.7.1.153) (Phosphatidylinositol 4,5-
--------------------------	--

Product Name: Recombinant Human PIK3CG Protein
Catalog #: PEH10592



bisphosphate 3-kinase 110 kDa catalytic subunit gamma) (PtdIns-3-kinase subunit p110-gamma) (p110gamma) (Phosphoinositide-3-kinase catalytic gamma polypeptide) (Serine/threonine protein kinase PIK3CG) (EC 2.7.11.1) (p120-PI3K)

Note

For Research Use Only , Not for Diagnostic Use