


Product Name: Recombinant Human UFM1 Protein
Catalog #: PEH10594



Summary

Name	UFM1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-GST
Accession #	P61960
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSKVSFKITL TSDPRLPYKV LSVPESTPFT AVLKFAAEF KVPAATSAII TNDGIGINPA QTAGNVFLKH GSELRIIPRD RVG

image



Background

Alternative Names Ubiquitin-fold modifier 1

Note

For Research Use Only , Not for Diagnostic Use