


Product Name: Recombinant Human PUM2 Protein
Catalog #: PEH10595



Summary

Name	PUM2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TB72
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IMVQRRSGQG FHGNSEVNAI LSPRSESGGL GVSMVEYVLS SSPADKLSR FRKGNFGTRD AETDGPEKGD QKGKASPFEE DQNRDLKQGD DDDSKINGRG LPNGMDADCK DFNRTPGSRQ ASPTEVVERL GPNTNPSEGL GPLNPNTANK PLVEEFSNPE TQNLDAEQV GLESLQFDYP GNQVPMDSG ATVGLFDYNS QQQLFQRTNA

image



Background

Alternative Names Pumilio homolog 2 (Pumilio-2)

Note

Product Name: Recombinant Human PUM2 Protein
Catalog #: PEH10595



For Research Use Only , Not for Diagnostic Use