


Product Name: Recombinant Human ZMIZ1 Protein
Catalog #: PEH10596



Summary

Name	ZMIZ1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9ULJ6
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PGPRPQSLPI QNIKRYPGE PNYGNQQYGP NSQFPTQPGQ YPAPNPPRPL TSPNYPGQRM PSQPSSGQYP PPTVNMGQYY KPEQFNGQNN TFSGSSYSNY SQGNVNRPPR PVPVANYPHS PVPGNPTPPM TPGSSIPPYL SPSQDVKPPF

image



Background

Alternative Names	Zinc finger MIZ domain-containing protein 1 (PIAS-like protein Zimp10) (Retinoic acid-induced protein 17)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use