


Product Name: Recombinant Human PAK6 Protein
Catalog #: PEH10605



Summary

Name	PAK6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NQU5
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TRESSLKRRRL FRSMFLSTAA TAPPSSSKPG PPPQSKPNSS FRPPQKDNPP SLVAKAQLSLP SDQPVGTFSP LTTSDTSSPQ

image



Background

Alternative Names	Serine/threonine-protein kinase PAK 6 (EC 2.7.11.1) (PAK-5) (p21-activated kinase 6) (PAK-6)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use