


Product Name: Recombinant Human SAMD4A Protein
Catalog #: PEH10608



Summary

Name	SMAG1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UPU9
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EDRTSTSFEGG QNRGRSDSVD YGQTHYYHQR QNSDDKLNGW QNSRDSGICI NASNWQDKSM GCENGHVPLY SSSSVPTTIN TIGTSTSTIL SGQAHHSPLK RSVSLTPPMN

image



Background

Alternative Names	Protein Smaug homolog 1 (Smaug 1) (hSmaug1) (Sterile alpha motif domain-containing protein 4A) (SAM domain-containing protein 4A)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use