


Product Name: Recombinant Human PKP4 Protein
Catalog #: PEH10615



Summary

Name	PKP4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q99569
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SKPRVSDAVQ PNNYLIRTEP EQGTLYSPEQ TSLHESEGL GNSRSSTQMN SYSDSGYQEA GSFHNSQNVK KADNRQQHSF IGSTNNHVVR NSRAEGQTLV QPSVANRAMR RVSSVPSRAQ SPSYVISTGV SPSRGLRTS LGSGFGSPSV TDPRLNPSA YSSTTLPAAR AASPYSQRPA SPTAIRRIGS VTSRQTSNPN GPTPQQYQTTA RVGSPLTLTD AQTRVASPSQ GQVGSSSPKR

image



Background

Alternative Names Plakophilin-4 (p0071)

Note

Product Name: Recombinant Human PKP4 Protein
Catalog #: PEH10615



For Research Use Only , Not for Diagnostic Use