


Product Name: Recombinant Human UBASH3B Protein
Catalog #: PEH10616



Summary

Name	UBS3B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	Q8TF42
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DFAAEAASKT EVHVEPHKKQ LHVTLAYHFQ ASHLPTLEKL AQNIDVKLGC DWVATIFSRD IRFANHETLQ VIYPYTPQND DELELVPGDF IFMSPMEQTS TSEGWIYGTS LTTGCSGLLP ENYITKADEC STWIFHGSYS ILNTSSSNSL TFGDGVLERR PYEDQGLGET TPLTIICQPM QPLRVNSQPG PQKRCLFVCR HGERMDVVFG KYWLSQCFDA KGRYIRTNLN MPHSLPQRSG GFRDYEKDAF

image



Background

Alternative Names	Ubiquitin-associated and SH3 domain-containing protein B (EC 3.1.3.48) (Cbl-interacting protein p70) (Suppressor of T-cell receptor signaling 1) (STS-1) (T-cell ubiquitin ligand 2) (TULA-2) (Tyrosine-protein phosphatase STS1/TULA2)
--------------------------	---

Product Name: Recombinant Human UBASH3B Protein
Catalog #: PEH10616



Note

For Research Use Only , Not for Diagnostic Use