


Product Name: Recombinant Human GAPDHS Protein
Catalog #: PEH10620



Summary

Name	G3PT_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	O14556
Host	E.coli
Species	Human
Molecular Mass	77 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NGFGRIGRLV LRACMEKGVK VVAVNDPFID PEYMVYMFKY DSTHGRYKGS VEFRNGQLVV DNHEISVYQC KEPKQIPWRA VGSPYVVEST GYVLSIQAAS DHISAGAQRV VISAPSPDAP MFVMGVNEND YNPGSMNIVS NASCTTNCLA PLAKVIHERF GIVEGLMTTV HSYTATQKTV DGPSRKAWRD GRGAHQNIIP ASTGAAKAVT KVIPELKGKL TGMAFRVPTP DVSVVDLTCR LAQPAPYSAI KEAVKAAAKG PMAGILAYTE DEVVSTDFLG DTHSSIFDAK AGIALNDNFV KLISWYDNEY

image



Background

Alternative Names	Glyceraldehyde-3-phosphate dehydrogenase, testis-specific (EC 1.2.1.12) (Spermatogenic cell-specific glyceraldehyde 3-phosphate dehydrogenase 2)
--------------------------	---

Product Name: Recombinant Human GAPDHS Protein
Catalog #: PEH10620



(GAPDH-2) (Spermatogenic glyceraldehyde-3-phosphate dehydrogenase)

Note

For Research Use Only , Not for Diagnostic Use