


Product Name: Recombinant Human GAPDHS Protein
Catalog #: PEH10621



Summary

Name	G3PT_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	O14556
Host	E.coli
Species	Human
Molecular Mass	78 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SVARELTVGI NGFGRIGRLV LRACMEKGVK VVAVNDPFID PEYMVYMFKY DSTHGRYKGS VEFRNGQLVV DNHEISVYQC KEPKQIPWRA VGSPYVVEST GVYLSIQAAAS DHISAGAQRV VISAPSPDAP MFVMGVNEND YNPGSMNIVS NASCTTNCLA PLAKVIHERF GIVEGLMTTV HSYTATQKTV DGPSRKAWRD GRGAHQNIIP ASTGAAKAVT KVIPELKGKL TGMAFRVPTP DVSVDLTCR LAQPAPYSAI KEAVKAAAKG PMAGILAYTE DEVVSTDFLG DTHSSIFDAK AGIALNDNFV KLISWYDNEY

image



Background

Alternative Names	Glyceraldehyde-3-phosphate dehydrogenase, testis-specific (EC 1.2.1.12) (Spermatogenic cell-specific glyceraldehyde 3-phosphate dehydrogenase 2)
--------------------------	---

Product Name: Recombinant Human GAPDHS Protein
Catalog #: PEH10621



(GAPDH-2) (Spermatogenic glyceraldehyde-3-phosphate dehydrogenase)

Note

For Research Use Only , Not for Diagnostic Use