


Product Name: Recombinant Human PTF1A Protein
Catalog #: PEH10626



Summary

Name	PTF1A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7RTS3
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ELLADEQAEV EFLSHQLHEY CYRDGACLLL QPAPPAAPLA LAPPSSGGLG EPDDGGGGGY CCETGAPPGG FPYSPGSPPS CLAYPCAGAA VLSPGARLRG LSGAAAAAAR RRRRVSEAE LQQLRQAANV RERRRMQSIN DAFEGLRSHI PTLPYEKRLS KVDTLRLAIG YINFLSELVQ ADLPLRGGGA GCGGPGGGG RLGGDSPGSQ

image



Background

Alternative Names	Pancreas transcription factor 1 subunit alpha (Class A basic helix-loop-helix protein 29) (bHLHa29) (Pancreas-specific transcription factor 1a) (bHLH transcription factor p48) (p48 DNA-binding subunit of transcription factor PTF1) (PTF1-p48)
--------------------------	---

Product Name: Recombinant Human PTF1A Protein
Catalog #: PEH10626



Note

For Research Use Only , Not for Diagnostic Use