


Product Name: Recombinant Human GULP1 Protein
Catalog #: PEH10628



Summary

Name	GULP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UBP9
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KTWMHTPEAL SKHFIPYNAK FLGSTEVEQP KGTEVVRDAV RKLKFARHIK KSEGQKIPKV ELQISYGVK ILEPKTKEVQ HNCQLHRISF CADDKTDKRI FTFICKDSES NKHLCYVFDS EKCAEEITLT IGQAFDLAYR KFLESGGKDV ETRKQIAGLQ KRIQDLETEN MELKNKVQDL ENQLRITQVS APPAGSMTPK

image



Background

Alternative Names	PTB domain-containing engulfment adapter protein 1 (Cell death protein 6 homolog) (PTB domain adapter protein CED-6) (Protein GULP)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human GULP1 Protein
Catalog #: PEH10628

