


Product Name: Recombinant Human TGM3 Protein
Catalog #: PEH10631



Summary

Name	TGM3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q08188
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LYDNTTGKQW KNSVNSHTIG RYISTKAVGS NARMDVTDKY KYPEGSDQER QVFQKALGKL KPNTPFPAATS SMGLETEEQE

image



Background

Alternative Names	Protein-glutamine gamma-glutamyltransferase E (EC 2.3.2.13) (Transglutaminase E) (TG(E)) (TGE) (TGase E) (Transglutaminase-3) (TGase-3) [Cleaved into: Protein-glutamine gamma-glutamyltransferase E 50 kDa catalytic chain; Protein-glutamine gamma-glutamyltransferase E 27 kDa non-catalytic chain]
--------------------------	--

Note

Product Name: Recombinant Human TGM3 Protein
Catalog #: PEH10631



For Research Use Only , Not for Diagnostic Use