


Product Name: Recombinant Human LPIN2 Protein
Catalog #: PEH10632



Summary

Name	LPIN2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92539
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GINQATLSGC IDVIVVQQQD GSYQCSPFHV RFGKLGVLRS KEKVIDIEIN GSAVDLHMKL GDNGEAFVVE ETEEEYEKLP AYLATSPIPT EDQFFKDIDT PLVKSGGDET PSQSSDISHV LETETIFTPS SVKKKKRRRK KYKQDSKKEE QAASAAAEDT CDVGVSDDDD KGAQAARGSS NASLKEEECK EPLLFHSGDH

image



Background

Alternative Names Phosphatidate phosphatase LPIN2 (EC 3.1.3.4) (Lipin-2)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human LPIN2 Protein
Catalog #: PEH10632

