

**Product Name: Recombinant Human DYRK3 Protein**  
**Catalog #: PEH10633**



## Summary

<b>Name</b>	DYRK3_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	O43781
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	13 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	TKCPPCSNVL CNPSEPPPPR RLNMTTEQFT GDHTQHFLDG GEMKVEQLFQ EFGNRKSNTI QSDGISDSEK CSPTVSQGKS SDCLNTVKS N SSSKAPKVVP LTPEQALKQY

## image



## Background

**Alternative Names** Dual specificity tyrosine-phosphorylation-regulated kinase 3 (EC 2.7.12.1)  
(Regulatory erythroid kinase) (REDK)

## Note

For Research Use Only , Not for Diagnostic Use