


Product Name: Recombinant Human PIK3C2B Protein
Catalog #: PEH10635



Summary

Name	P3C2B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O00750
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QNADPSLISW DEPGVDFYSK PAGRRTDLKL LRGLSGSDPT LNYNSLSPQE GPPNHSTSQG PQPGSDPWPK

image



Background

Alternative Names	Phosphatidylinositol 4-phosphate 3-kinase C2 domain-containing subunit beta (PI3K-C2-beta) (PtdIns-3-kinase C2 subunit beta) (EC 2.7.1.154) (C2-PI3K) (Phosphoinositide 3-kinase-C2-beta)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use